

Development of Sum Frequency Generation Imaging Microscopy


Dr. Zenonas Kuprionis
EKSPLA UAB

Outline

- ❑ **Basic principles of SFG spectroscopy**
- ❑ **Features of SFG**
- ❑ **SFG spectrometer**
- ❑ **Challenges of SFG imaging**
- ❑ **Possible solutions**
- ❑ **EKSPILA SFG microscope**

Basic principles of Vibrational SFG Spectroscopy

- Sum Frequency Generation (SFG) is a second-order nonlinear optical process, which is used to generate **a vibrational spectrum of molecules at an interface**.
- SFG is generated through overlapping of infrared (ω_{IR}) and visible (ω_{VIS}) laser beams to produce an output at the sum frequency ($\omega_{SFG} = \omega_{IR} + \omega_{VIS}$).
- SFG allowed in media without inversion symmetry. At surfaces and interfaces inversion symmetry is necessary broken, so SFG can be used as an effective surface probe.
- Different polarization combinations of pump beams and analysis of SFG beam polarization provides information about orientation of adsorbate vibrational modes.



$$I(\omega_{SFG}) \sim |\chi^{(2)}|^2 I(\omega_{VIS}) I(\omega_{IR})$$

$$\chi^{(2)} = \chi^{(2)}_{NR} + \chi^{(2)}_R$$

$$\chi^{(2)}_R \sim A_{\nu} M_{\nu,g} / (\omega_{IR} - \omega_0 - i\gamma)$$

SFG as an Analytical Technique

Method	Surface/Interface Monolayer Specificity	Good Spatial and Spectral Resolution	Structural Information	Functional Group Orientation	Suitable for Liquid/Solid, Gas/Solid, Liquid/Gas, Liquid/Liquid Interfaces	in situ Measurement	Effective at Polymer-Liquid Interface	Non-destructive
SFG	●	●	●	●	●	●	●	●
Reflection Infrared Spectroscopy		●	●	●*	●	●*		●
Attenuated Total Reflection Infrared Spectroscopy		●	●	●*	●	●*	●*	●
Raman Spectroscopy		●	●	●*		●*		
Contact Angle Measurement							●	●
Neutron Reflection	●	●	●					
X-Ray Photoelectron Spectroscopy (XPS)	●	●	●					
Surface Light Scattering	●	●				●	●	●
Secondary Ion Mass Spectroscopy (SIMS)	●	●	●					

* Limited conditions apply


Table 1. Comparison of SFG to other conventional surface analysis techniques

(Polymer Technology Group, Inc.)

Powerful tool for *in-situ* investigation of surfaces and interfaces:

- Highly surface-specific
- Sensitive to sub-monolayers
- Applicable to all interfaces accessible by light
- Molecular selectivity
- Remote sensing in hostile environment
- Surface dynamics investigation


Vibrational Sum Frequency Generation Spectrometer



- Picosecond Nd:YAG laser
- Harmonics generation unit
- Optical parametric generator
- Monochromator
- Signal detectors
- Data acquisition system
- Controlling software
- Beam delivery optics
- Guiding beam for system alignment

SFG imaging advantages

- High interface specificity
- Monolayer sensitivity
- Chemical contrast
- Structural information
- Spatial distribution
- Answering questions of sample inhomogeneity



Difficulties:


- Low efficiency
- Sensitivity of samples to damage
- Oblique incidence of pump beams

Possible solutions of SFG imaging

- Raster scanning
- Near-field microscopy
- Far-field imaging

Near-field microscopy

- High resolution (200 nm)
- Low coupling efficiency to fiber
- No results with vibrational resonances



Richard D. Schaller and Richard J. Saykally
Department of Chemistry, UC, Berkeley

Samples of near-field microscopy

Langmuir 2001, 17, 2055–2058


Near-Field Infrared Sum-Frequency Generation Imaging of Chemical Vapor Deposited Zinc Selenide

Richard D. Schaller and Richard J. Saykally*

Department of Chemistry, University of California, Berkeley, California 94720-1460

SFG NSNOM images of MCD ZnSe

White features are interpreted as strain patterns within ZnSe crystal lattice.



Far-field imaging


Problem of SFG geometry:

Optimal pump beams incidence angles 50 – 60 deg

Undistorted oblique imaging?

Solution:


Detection via an intermediate image on the surface of a blazed grating



EKSPLA SFG imaging microscope (joint project with Steven Baldelli group, HU)




Specifications

Resolution	<5 μm
Field of view	300 μm
Magnification	20, 10
Wavelength range	2000 – 3800 cm^{-1}



An image of ODT covered etched gold sample. The first lines are spaced, in order left -right 20, 5, 1, 2, and 8 micrometers. Image taken at $\omega_{\text{IR}}=2875 \text{ cm}^{-1}$ and acquired for 5000 laser shots.



Origin of chemical contrast



2875 cm⁻¹ - represent resonant CH₃

2955 cm⁻¹ - represent non-resonant

3060 cm⁻¹ - represent resonant phenyl



hexadecane-1-thiol

16-phenylhexadecane-1-thiol

Conclusions

- ❑ **SFG imaging microscope with few μm resolution was developed and tested**
- ❑ **Microscope was successfully used for imaging of:**
 - ✓ Microcontact-printed mixed self-assembled monolayers
 - ✓ CO deposition on PT
 - ✓ Polymer blends